


Facile & Reliable Screening

Screening Audiometer Ultimate



Pejvak Ava Corp.
Audiology Equipment



Ultimate Screening Audiometer

Ultimate Technical Specifications

Standards :

Tone audiometer type:4
Calibration:ISO389

Channels : One channel.
3 Mode function.

Frequencies and Maximum Hearing Levels:

Hz \ Mode	1	2	3
125	*	*	*
250	80	*	*
500	100	100	50
750	100	*	*
1000	100	100	50
1500	100	*	*
2000	100	100	50
3000	100	*	*
4000	100	100	50
6000	80	*	*
8000	80	*	*
Min HL	-10	-10	0 dB
dB Steps	5	5	10 dB

Input/ Output: Tone/ Left& Right.

Presentations: Manual,
single pulses.

Frequency Range: 250-8000 Hz.

Frequency Resolution:1/2 octave steps.

Modulation:Warble:5Hz

Attenuator:
Totally Click free, -10 to 100 dB HL in 5 &10dB Steps.

Tone Switches:
Silent touch tone switch.

Patient Response:
One patient response button,

Reset Mode:
For Quick audiometry.

Test Types:
Tone: Manual, single pulse,
2 Mode Automatic Screening Audiometry.

Free Field:
2x12W external amplifier and two speakers.
(Optional)

Communication:
No Available.

Dimensions: LxWxH: 29x24.5x5 cm.

Weight: 2000 gr.

Power: AC 50-60 Hz.180-240 V
Auto Power Off System (STD)
by alarm mode

Consumption: Max. 5 VA

Included Parts:

- Bilsom Screening Audiometric Headset
(Audiocup Noise Excluders)
- 1 Patient Response Button
- Power Cable 240V
- Operation Manual
- Soft Carrying case

Options:

- Ap12 Built-in 2x12 watt Power Amplifier
for Free Field
- AP12 Loudspeaker
- Hard Carrying case

